LABOSHAKE

LS 2/5, RO 2/5

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Please read this instruction manual with care before you start operating the system!

Please observe the safety instructions of this manual, marked with  in order to avoid any dangerous handling.

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Contents

	Safety instructions	16
1.	General remarks	17
	Warranty conditions	17
	Transport damage	17
2.	Assembly	18
	Delivery	18
	Electrical connection	18
	Flasks	18
	Technical data L A B O S H A K E	18
3	Accessories/Shaking attachments	18
	Universal attachments	19
	Fig. 1: Universal attachment UA	19
	Carrier plates	19
	Fig. 2: Carrier plates and two-storey frame	21
	Fig. 3: Holder for separating funnels	21
4	Setting into operation	22
	Keyboard	22
	Fig. 4: Keyboard	22
5	Programming	23
	Set up speed	23
	Select time-mode	23
	Set up time	23
6	Start	24
7	RS485-Interface	24
	Synchronous mode	25
8	Maintenance	25
9	Spare parts	25

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Safety instructions

Make sure that no liquid gets into contact with cable connections or the interior of the electrical parts of the unit!
Danger of electric shock!

For repairs on electrical, electronic or mechanical parts please always contact the service engineer of your agent. Don't do repairs by yourself!

Always switch off the apparatus at the mains and pull out the plug before opening! Danger of electric shock!

When switched off please wait until the machine has stopped completely.
Risk of injury!

Take care when handling glass parts and follow the national safety regulations concerning the handling of glass parts.

In order to ensure stability it has to be placed on to an even, horizontal area. Please ensure that there is enough safety distance between other appliances.

The shaker must **not be run in damp or hazardous location**. The maximum humidity allowed is 80%, the maximum room temperature must not exceed 40° C.

1. General remarks

Warranty Conditions

The LABOSHAKE is manufactured following the high quality guidelines of DIN EN ISO 9001. On the basis of the C. GERHARDT conditions of warranty our products are guaranteed for 1 year, as long as the apparatus is used in accordance with the instructions mentioned in this manual.

Transport damage

LABOSHAKE shakers are delivered in purpose built packing. So transport damage is nearly impossible. In case of any damage, please notify the carrier (post, rail, road) immediately and try to obtain an expert's report.

2. Assembly

Delivery

The delivery of the basic model comprises the connection cable and 4 rods, 150 mm long, to fasten the different shaking attachments.

Electrical connection

Connect the appliance to a shockproof socket by the cable enclosed. Please check the nominal voltage on the type plate before the connection to the mains. Before you first take the Laboshake into operation make sure that the excess current switch at the rear is in position 1.

Flasks

All types of LABOSHAKE can be equipped with almost all kinds of flasks.

Technical data

LABOSHAKE LS 2/5, RO 2/5

	LS 2 / RO 2	LS 5 / RO 5
Nominal voltage	230 V AC / 50 Hz	230 V AC / 50 Hz
Power consumption	80 W	80 W
Amplitude	20 mm	50 mm
Microprocessor	Single-chip (MCU) mit EEPROM	
Speed frequency	20-350 min ⁻¹	20-200 min ⁻¹
Time program (minute mode)	0 to 999 min or infinity	0 to 999 min or infinity
Time program (hour mode)	0 to 99 h und 0 to 59 min or infinity	0 to 99 h und 0 to 59 min or infinity
Load capacity	15 kg	30 kg
Shaking platform	370x300 mm	540x370 mm
Dimensions (LxDxH)	410x380x120 mm	610x490x120 mm
Weight	17 kg	27 kg

3. Accessories/Shaking attachments

The different models of shaking attachments are delivered as kit sets, easy to assemble. They are fixed on the 4 rods of the basic models.

Universal

attachment

These multi-purpose attachments are supplied as a kit set, easy to assemble. They consist of two fastening bars made of stainless steel, which are fixed by quick-clamps horizontally at any height onto the 4 rods to the basic models.

To these fastening bars the clamping rolls, covered with Neoprene-cushioning, are mounted for the safe clamping of a wide range of vessels in the upright position.

Universal attachment UA 2 for small shakers LS/RO 2

kit set, easy to assemble, consisting of

2 fastening bars 375 mm UB 25, No. 4213

2 clamping rolls 305 mm US 20, No. 4215

Universal attachment UA 5 for large shakers LS/RO 5

kit set, easy to assemble, consisting of


2 fastening bars 545 mm UB 50, No. 4214

3 clamping rolls 375 mm US 25, No. 4216

For special shaking requirements additional fastening bars and clamping rolls may be fixed on the attachment.

For shaking large bottles in horizontal position 2 additional fastening bars (No. 4212 or No. 4213) and clamping rolls (No. 4216 or 4217) can be mounted. With these additional parts nearly all kinds of shaking tasks can be performed

Fig. 1: Universal attachment UA with additional parts

Carrier plates

Carrier plates are available as universal plates for free assembly with clip racks for Erlenmeyer flasks of choice and finished plates taking a fixed quantity of Erlenmeyer flasks.

Universal plate TB 20 for small shakers LS/RO 2

330x300 mm, stainless steel, with a grid of 5 mm holes spaced at 20 mm intervals (18x15 = 270 holes). It is mounted directly onto the rods and fixed automatically.

Universal plate TB 50 for large shakers LS/RO 5

500x380 mm, stainless steel with a grid of 5 mm holes spaced at 20 mm intervals (25x19 = 475 holes). It is supplied with 2 guide rails, which are screwed on the rods. The carrier plate slides on these guide rails and is fixed by lateral screws.

To be equipped with clip racks (No. 4233 - 4239) of your choice.

Clip racks

4233	TK	25	Erlenmeyer	25 ml	TB 20:	28 ea.	TB 50:	54 ea.
4234	TK	50	Erlenmeyer	50 ml	TB 20:	28 ea.	TB 50:	48 ea.
4235	TK	100	Erlenmeyer	100 ml	TB 20:	16 ea.	TB 50:	33 ea.
4236	TK	250	Erlenmeyer	250 ml	TB 20:	9 ea.	TB 50:	20 ea.
4237	TK	500	Erlenmeyer	500 ml	TB 20:	8 ea.	TB 50:	12 ea.
4238	TK	1000	Erlenmeyer	1000 ml	TB 20:	5 ea.	TB 50:	8 ea.
4239	TK	2000	Erlenmeyer	2000 ml	TB 20:	2 ea.	TB 50:	5 ea.

Finished plate TB 21 for small shakers

similar to TB 20, but equipped with clamping spring assembly (V2A) taking 16 Erlenmeyer 50 up to 100 ml. It is mounted directly on the rods and fixed automatically.

Finished plate TB 51 for large shakers

similar to TB 50, but equipped with clamping spring assembly (V2A) taking 30 Erlenmeyer 50 up to 100 ml. It is supplied with 2 guide rails which are screwed on the rods. The carrier plate slides on these guide rails and is fixed by lateral screws.

Finished plate TB 52 for large shakers


similar to TB 51, but for taking 20 Erlenmeyer 200-300 ml.

Three-storey frame for large shakers

for shaking 3 carrier plates at a time. Guide rails take in 3 levels of universal plates TB 50 or finished carrier plates TB 51/52 of your choice.

Consisting of 4 rods EE 50, no. 40690, 460 mm replacing the shorter rods of the basic shaker and 8 spacing bars ED 50, no. 40662, 190 mm. Supplied as a kit set, easy to assemble.


Fig. 2: Carrier plates and two-storey frame

Two-piece holder HS 25, no. 4250 made of V2A to be screwed on universal plates TB 20 or TB 50. The maximum loading depends on the size of the separating funnel:

TB 20:	max.3 separating funnels up to 500 ml
TB 50:	max.2 separating funnels up to 2000 ml
	max.3 separating funnels up to 1000 ml
	max.4 separating funnels up to 500 ml
	(with three-storey frame EA 3 max. size 250 ml)

Fig. 3: Holder for separating funnels

4. Setting into operation

Excess current switch


Before you first take the Laboshake into operation make sure that the excess current switch at the rear is in position 1.

Keyboard

All C. Gerhardt shakers are equipped with modern electronic devices (single-chip MCU) that allow a precise adjustment of speed and time by the keyboard.

To extend durability of the display elements, they will be dimmed a few seconds after last key is pressed.

Fig. 4: Keyboard

- | | | |
|---|---------------------------------------|--|
| 1 | Prog - key/
-display | Key / display of speed |
| 2 | Time - key/
-display | Key / display of time |
| 3 | + | Increase speed and time |
| 4 | - | Decrease speed and time |
| 5 | Run | Soft start of the shaker:
After pressing the key „Run“, the machine starts running and increases speed softly until set speed is reached.
Quick start of the machine:
Running the key „Run“ and keeping it pressed forces the machine to increase speed rapidly until the set speed is reached. |
| 6 | Stop | Stopping the machine
After the key „Stop“ is pressed the shaker decreases the speed to zero. The selected values of speed and time remain stored when the machine is restarted with the key „Run“. |
| 7 | | Mains switch |

5. Programming

Set up speed

The PROG display shows the actual value when the power is turned on. Pressing „Prog“ the display shows the set value, indicated by a flashing dot. It can be changed with the keys „+“ or „-“.

Select time - mode (hour or minute)

There are two modes of displaying the time:

- Minute-mode: Ranges from 0-999 min or infinity.
- Hour-mode: Ranges from 0-99 h and 0-59 min or infinity.

The default setting is minute-mode.

There are two possible ways to change from minute- to hour-mode:

1. **Temporary change:** Pressing the keys „Time“ and „Stop“ simultaneously changes from minute-mode to hour-mode and vice versa.

The set mode remains, until the shaker is turned off. When the shaker is turned on again it runs in default mode.

2. **Long-term change:** While turning on the shaker press the key „Time“ and release it after the beep-signal (if it is not the desired mode, repeat the last step).

The selected mode remains when the shaker is turned off and on again (it is still possible to change modes temporarily).

Set up time

Press the key „Time“.

The TIME-display changes anytime the key is pressed as follows:

(Minute-mode)

Minute-set value → minute-actual value.

(Hour-mode)

Minute-set value → minute-actual value → hour-set value → hour-actual value.

The set values are indicated by a flashing dot. They can be changed by pressing the keys „+“ or „-“.

If you are in the minute-mode, you can select from zero to 999 min or infinity. In the hour-mode you can select from 0 to 99 h and 0 to 59 min or infinity. If you select infinity, the display shows „---“. This can be achieved by increasing the time to more than 999 or decreasing it below 000. If you change the set time during operation, the elapsed time is subtracted from the new set time.

Example: Set time was 50 min, the shaker has been running for 5 min, the new time value has been set to 10 min. The shaker will have a remaining time of 5 min. If the old set time was infinity, the elapsed time is calculated in the same way.

If you set the time in hour mode > 16h 39min (> 999min), the display will show „]]“ in minute-mode. This is due to 3-digit displays, which will over-flow at time values > 999. In this case any changes in time values can be made only in the hour mode or if the shaker is stopped. The display „]]“ turns to normal when the elapsed time has decreased to 999 min.

6. Start

CAUTION: Wait until the machine has stopped completely before you load or unload it.
Risk of injury!

To keep up the durability of the mechanical parts, never exceed the maximum load capacity (see technical data LABOSHAKE).

Start

The machine is turned on by pushing the mains switch. Then it is started with the key „Run“, provided the time value is > 0 min. Values of speed and time from earlier programs are restored and used again. They can be changed anytime. The soft start protects from heavy motion due to false programming. You can interrupt the shaker operation any time by pressing the key „Stop“ and correct the set values if necessary.

7. The RS485-Interface

The interface of the shaker allows you:

- a) to handle up to 32 devices synchronously (synchro-mode), **preferably of the same kind**. This means, one shaker is programmed as usual and all other shakers connected to the bus-system are programmed with this shaker data. The corresponding number of cables, order no. 11243, and one termination plug, order no. 4265, must be used to provide a connection between the devices.
- b) to operate several devices of **different** kinds using a PC and the Gerbus[®]-applications. This enables you to gain full control of all devices involved. The Gerbus[®]-applications can be received on request.

Synchro-mode

To operate in synchro-mode you need at least two machines **(preferably of the same type)** which have to be connected via RS485-cable. The terminator is plugged to one of the machines.

Then turn on all machines involved.

The interface on all machines are activated by pressing the keys „Time“ and „+“ simultaneously. The time display shows a „□“ in the lower part of the TIME-display to indicate the machine is „ON-Line“.

Repeat the last step for all machines involved.

Now you select one of the machines to be the „Master“, the one that gets all the key commands. This is done by pressing „Prog“ and „+“ simultaneously. The „Master“ is indicated by a „□“ in the upper part of the TIME-display.

Now all the key commands and values you set on the „Master“, are received and copied by all machines connected to the bus-system.

The synchro-mode can be used in minute- and hour-mode.

You deactivate the „Master“ by pressing the keys „Prog“ and „-“ simultaneously.

To go „Off-Line“ completely, press the keys „Time“ and „-“ or turn the machine off and on again.

This goes for all machines connected in the bus-system.

8. Maintenance

The shakers are almost maintenance-free. Only the driving belts and possibly the self-adhesive feet wear, due to rubbing, and need checking. The self-adhesive feet can be easily replaced by yourself. The driving belts have to be replaced by authorized service personnel. Pull out the mains plug and remove the shaking platform by lifting the rubber matting and loosening the screws.

If there are problems with the electronic control system, please contact the service personnel of your local agent.

9. Spare parts

Cat.-no.	Description	LS 2	LS 5	RO 2	RO 5
CASING					
15902	Self-adhesive foot, M8	4	4	4	4
11188	Keyboard SM, long		1		1
11189	Keyboard SM, short	1		1	
11750	Main switch, green	1	1	1	1
15370	Excess current switch	1	1	1	1
1311	Connection cable	1	1	1	1
14401	Power socket	1	1	1	1
14626	Drive-belt, D4x672 mm	1		1	
14621	Drive-belt, D4x360 mm		2		2
14622	Drive-belt, D4x820 mm		2		2

Cat.-no.	Description	LS 2	LS 5	RO 2	RO 5
SHAKING PLATFORM					
40621	Driver 100x70x10 m	1	1		
40622	Driver D75x10			1	1
40641	Running surface coat 80x10x2,5 mm	1	1		
12402	Rubber ring black, D32/6x3 m	4	4	4	4
40653	Corner bolt D10x155 m M6	4	4	4	4
MOTOR UNIT					
10501	Flex-lock rubber piece M6	4	4	4	4
11197	Optoelectronic sensor with plug	1	1	1	1
12220	External rotor motor M2E068-CF0101	1	1	1	1
12225	Motor condensator MFB-MP2/400	1	1	1	1
40616	Motor disc D40	1	1	1	1
40657	Hole disc D111	1	1	1	1
FLYWHEEL BEARING					
15604	Ball bearing SKF, 6201-2Z/QE6	1	1	1	1
15614	Ball bearing SKF, 6203-2ZNR/QE6	2	2	2	2
BEARING BLOCK 30x30x78					
15660	Track roller KLR-U 08.28.11.2Z	1	1		
40642	Slide way cover 50x10x2,5	2			
40643	Slide way cover 80x10x2,5		2		
BOARD FIXING DEVICE					
11184	SM-Flat strip 10/10 pol.	1	1	1	1
11195	SM-Power PCB 220 V	1	1	1	1
11198	SM-Power PCB 240 V	1	1	1	1
11221	SM-Control PCB	1	1	1	1
INTERMEDIATE DISC					
15614	Ball bearing SKF,6203-2ZNR/QE6		2		2